

ABX Minocal



LOT

CX 489

CAL



(Exp.) 2024-01-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500 V1.0 to V1.4	H1500 H2500 Since V1.7							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.55	8.50							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.53							± 0.06
HB HGB	g/dl	13.3	13.3							± 0.2
	g/l	133	133							± 2
	mmol/l	8.26	8.26							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	244	236							± 10
VGM MCV	µm ³ ; fl	83.0	84.2							± 4.0

1300116192
Ref: TEMP-0387 Rev.50