

LOT		CX 487		CAL		(Exp.) 2023-11-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN								TOLERANCES TOLERANCE
		H1500	H1500								
		H2500	H2500								
		V1.0 to V1.4	Since V1.7								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.88								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57								± 0.06
HB HGB	g/dl	13.3	13.2								± 0.2
	g/l	133	132								± 2
	mmol/l	8.26	8.20								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	245	247								± 10
VGM MCV	μm ³ ; fl	82.0	82.0								± 4.0

1300116192

Ref. TEMP-0387 Rev.50

 1300116192
 Ref. TEMP-0387 Rev.50