

ABX Minocal



LOT

CX 484

CAL



(Exp.) 2023-08-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56								± 0.06
HB HGB	g/dl	13.6								± 0.2
	g/l	136								± 2
	mmol/l	8.45								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	254								± 10
VMP MPV	µm ³ ; fl	N/A								N/A
VGM MCV	µm ³ ; fl	84.0								± 4.0

1300060884
Ref: TEMP-0387 Rev.49