

# ABX Minocal



**LOT**

**CX 481**

**CAL**



**(Exp.) 2023-05-05**  
(YYYY-MM-DD)

Rev **2**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.55</b>								± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.41</b>								± 0.06
HB HGB	g/dl	<b>12.9</b>								± 0.2
	g/l	<b>129</b>								± 2
	mmol/l	<b>8.01</b>								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>254</b>								± 10
VMP MPV	µm <sup>3</sup> ; fl	N/A								N/A
VGM MCV	µm <sup>3</sup> ; fl	<b>82.5</b>								± 4.0

1300060884  
Ref: TEMP-0387 Rev.48