

ABX Minocal



LOT CX 475
Rev 1

CAL

 (Exp.) 2022-11-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	9.00								± 0.20
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.58								± 0.06
HB HGB	g/dl	13.6								± 0.2
	g/l	136								± 2
	mmol/l	8.45								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLA PLT	10 ³ /mm ³ ; 10 ⁹ /l	244								± 10
VMP MPV	µm ³ ; fl	N/A								N/A

1300060884
L-0387 Rev.48

ABX Minocal



VGM	MCV	μm^3 ; fl	83.0																		± 4.0
-----	-----	----------------------	-------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	-----------

Ref: TEMP