

ABX Minocal

LOT		CX 488		CAL		 (Exp.) 2023-12-05 (YYYY - MM - DD)						
Rev 1		ABX Lysebio										TOLERANCES TOLERANCE
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN		YUMIZEN								
		H1500 H2500	H1500 H2500									
		V1.0 to V1.4	Since V1.7									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.70	8.79								± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.43								± 0.06	
HB HGB	g/dl	13.1	13.1								± 0.2	
	g/l	131	131								± 2	
	mmol/l	8.14	8.14								± 0.12	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	236								± 10	
VGM MCV	μm ³ ; fl	82.5	82.5								± 4.0	

13001.16192

Ref: TEMP-0387 Rev.50

 1300116192
 Ref. TEMP-0387 Rev.50