

ABX Minocal



LOT		CX 489		CAL		 (Exp.) 2024-01-05 (YYYY - MM - DD)							
		Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN		YUMIZEN								TOLERANCES TOLERANCE	
		H1500	H1500										
		H2500	H2500										
		V1.0 to V1.4	Since V1.7										
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.55	8.50									± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.53									± 0.06	
HB HGB	g/dl	13.3	13.3									± 0.2	
	g/l	133	133									± 2	
	mmol/l	8.26	8.26									± 0.12	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	244	236									± 10	
VGM MCV	µm ³ ; fl	83.0	84.2									± 4.0	

1300116192
Ref. TEMP-0387 Rev.50