

ABX Minocal



LOT		CX 472		CAL		(Exp.) 2022-08-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500									
		H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75									± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54									± 0.06
HB HGB	g/dl	13.5									± 0.2
	g/l	135									± 2
	mmol/l	8.38									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	245									± 10
VMP MPV	μm ³ ; fl	N/A									N/A
VGM MCV	μm ³ ; fl	85.0									± 4.0

 1300060884
 Ref: TEMP-0387 Rev.48