

LOT		CX 486		CAL		 (Exp.) 2023-10-05 (YYYY - MM - DD)					
Rev 2											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500	H1500								
		H2500	H2500								
		V1.0 to V1.4	Since V1.7								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.70	8.72							± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.50							± 0.06	
HB HGB	g/dl	13.2	13.2							± 0.2	
	g/l	132	132							± 2	
	mmol/l	8.20	8.20							± 0.12	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	241	241							± 10	
VGM MCV	μm ³ ; fl	82.5	83.0							± 4.0	

1300116192

Ref. TEMP-0387 Rev.50

 1300116192
 Ref. TEMP-0387 Rev.50