

ABX Minocal



LOT		CX 483		CAL		(Exp.) 2023-07-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55									± 0.06	
HB	HGB	g/dl	13.2									± 0.2	
		g/l	132									± 2	
		mmol/l	8.20									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	258									± 10	
VMP	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	81.0									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.49