

ABX Minocal

LOT		CX 476		CAL		(Exp.) 2022-12-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80									± 0.20		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55									± 0.06		
HB HGB	g/dl	13.3									± 0.2		
	g/l	133									± 2		
	mmol/l	8.26									± 0.12		
HT HCT	%	N/A									N/A		
	l/l	N/A									N/A		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	250									± 10		
VMP MPV	µm ³ ; fl	N/A									N/A		
VGM MCV	µm ³ ; fl	82.5									± 4.0		

 1300060884
 Ref: TEMP-0387 Rev.48