

ABX Minocal

LOT
CX 486

 Rev **2**
CAL

(Exp.) 2023-10-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500	H1500								
		H2500	H2500								
		V1.0 to V1.4	Since V1.7								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.70	8.72								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.50								± 0.06
HB HGB	g/dl	13.2	13.2								± 0.2
	g/l	132	132								± 2
	mmol/l	8.20	8.20								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	241	241								± 10
VGM MCV	μm ³ ; fl	82.5	83.0								± 4.0

 1300116192
 Ref. TEMP-0387 Rev.50