

LOT		CX 487		CAL		 (Exp.) 2023-11-05 (YYYY - MM - DD)						
Rev 1		ABX Lysebio										TOLERANCES TOLERANCE
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN									
		H1500	H1500									
		H2500	H2500									
		V1.0 to V1.4	Since V1.7									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.88								± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57								± 0.06	
HB HGB	g/dl	13.3	13.2								± 0.2	
	g/l	133	132								± 2	
	mmol/l	8.26	8.20								± 0.12	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	245	247								± 10	
VGM MCV	μm ³ ; fl	82.0	82.0								± 4.0	

1300116192

Ref: TEMP-0387 Rev.50

 1300116192
 Ref. TEMP-0387 Rev.50