

ABX Minotrol Retic



LOT RX 436
Rev 1

CONTROL

(Exp.) 2022-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.5			± 0.8	5.5			± 2.2	10.4			± 3.1
RETIC	%	1.82			± 0.91	5.47			± 1.37	12.43			± 2.49
RETIC	10 ⁶ /mm ³	0.064			± 0.032	0.193			± 0.048	0.434			± 0.087
	10 ⁹ /l	64			± 32	193			± 48	434			± 87
MFI MRI *	%	29.0			± 8.0	41.6			± 10.0	48.1			± 12.0
	-	0.290			± 0.080	0.416			± 0.100	0.481			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.51			± 0.30	3.52			± 0.30	3.49			± 0.30
PIC RANGE	-	13			± 8	13			± 8	13			± 8

* For research use only in Canada and USA

1300102875
Ref: TEMP-0382 Rev.42