

ABX Minotrol Retic



LOT RX 441
Rev 1

CONTROL

(Exp.) 2023-07-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.3			± 0.7	6.2			± 2.5	10.0			± 3.0
RETIC	%	1.60			± 0.80	5.27			± 1.32	10.71			± 2.14
RETIC	10 ⁶ /mm ³	0.056			± 0.028	0.182			± 0.046	0.367			± 0.073
	10 ⁹ /l	56			± 28	182			± 46	367			± 73
MFI MRI *	%	31.9			± 8.0	43.9			± 10.0	48.5			± 12.0
	-	0.319			± 0.080	0.439			± 0.100	0.485			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.49			± 0.30	3.46			± 0.30	3.43			± 0.30
PIC RANGE	-	19			± 8	19			± 8	19			± 8

1300102875
Ref: TEMP-0382 Rev.42

* For research use only in Canada and USA