

ABX Minotrol Retic



LOT RX 442
Rev 1

CONTROL

(Exp.) 2023-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500				H2500				H2500			
Méthode Manuelle / Manual RETIC	%	1.9			± 1.0	5.7			± 2.3	11.9			± 3.6
RETIC	%	1.43			± 0.72	5.22			± 1.31	11.51			± 2.30
RETIC	10 ⁶ /mm ³	0.049			± 0.025	0.180			± 0.045	0.398			± 0.080
	10 ⁹ /l	49			± 25	180			± 45	398			± 80
MFI MRI *	%	30.5			± 8.0	44.7			± 10.0	49.6			± 12.0
	-	0.305			± 0.080	0.447			± 0.100	0.496			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.45			± 0.30	3.44			± 0.30	3.46			± 0.30
PIC RANGE	-	19			± 8	19			± 8	19			± 8

1300102875
Ref: TEMP-0382 Rev.43

* For research use only in Canada and USA