

ABX Minotrol Retic



LOT RX 443
Rev 1

CONTROL

(Exp.) 2023-11-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500 V1.0 to V1.4				H2500 V1.0 to V1.4				H2500 V1.0 to V1.4			
Méthode Manuelle / Manual													
RET	%	1.4			± 0.7	5.5			± 2.3	11.8			± 3.3
RET	%	1.82			± 0.91	5.58			± 1.40	10.91			± 2.18
RET	10 ⁶ /mm ³	0.063			± 0.032	0.191			± 0.048	0.370			± 0.074
	10 ⁹ /l	63			± 32	191			± 48	370			± 74
MFI MRI *	%	32.3			± 8.0	42.1			± 10.0	46.0			± 12.0
	-	0.323			± 0.080	0.421			± 0.100	0.460			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.46			± 0.30	3.42			± 0.30	3.39			± 0.30
PIC RANGE	-	19			± 8	19			± 8	19			± 8

Ref: TEMP-0382 Rev.44 FRONT / RECTO 1300116193

* For research use only in Canada and USA