

ABX Minotrol Retic



LOT RX 445

CONTROL

(Exp.) **2024-03-05**
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500 HELO 2.0				H2500 HELO 2.0				H2500 HELO 2.0			
Méthode Manuelle / Manual													
RET	%	1.4			± 0.7	5.4			± 2.2	10.4			± 3.2
RET	%	1.50			± 0.75	5.03			± 1.26	9.69			± 1.94
RET	10 ⁶ /mm ³	0.052			± 0.026	0.175			± 0.044	0.333			± 0.067
	10 ⁹ /l	52			± 26	175			± 44	333			± 67
CRC	%	0.94			± 0.47	3.30			± 0.99	6.23			± 1.56
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.45			± 0.30	3.47			± 0.30	3.44			± 0.30
IRF	--	0.24			± 0.14	0.32			± 0.19	0.35			± 0.21
RHCC	pg	21.8			± 4.4	15.3			± 3.1	11.0			± 2.8

Ref: TEMP-0382 Rev.45 1300116193 BACK / VERSO