

LOT **BFI 0922**
Rev 1

CONTROL

 **(Exp.) 2022-11-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	BFTROL									
		CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE
		RANGE	RANGE	2	RANGE		RANGE	3			
YUMIZEN H1500	YUMIZEN H2500			YUMIZEN H1500	YUMIZEN H2500						
BFWBC	10 ⁶ /l	195	195			± 100	850	850			± 120
BFRBC	10 ⁶ /l	69460	69460			± 12000	476290	476290			± 50000

Ref: TEMP-5502 Rev.05 1300060764