

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

BFI 0523

Rev 1

CONTROL



(Exp.) 2023-07-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	BFTROL									
		CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE
		RANGE		RANGE			RANGE		RANGE		
YUMIZEN H1500	YUMIZEN H2500			YUMIZEN H1500	YUMIZEN H2500						
BFWBC	10 ⁶ /l	185	185			± 100	870	870			± 120
BFRBC	10 ⁶ /l	68000	68000			± 12000	465000	465000			± 50000

Ref: TEMP-5502 Rev.05 1300060764