

LOT BFI 2444
Rev 1

CONTROL

 (Exp.) 2024-01-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	BFTROL									
		CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE
		RANGE		RANGE			RANGE		RANGE		
		YUMIZEN H1500	YUMIZEN H2500	YUMIZEN H1500	YUMIZEN H2500		YUMIZEN H1500	YUMIZEN H2500	YUMIZEN H1500	YUMIZEN H2500	
		Since V1.7			Since V1.7						
BFWBC	10 ⁶ /l	237	237			± 95	908	908			± 227
BFRBC	10 ⁶ /l	77691	77691			± 19423	508674	508674			± 127169
BFMN	10 ⁶ /l	68	68			± 24	102	102			± 41
	%	28	28			± 10	11	11			± 4
BFPN	10 ⁶ /l	169	169			± 59	806	806			± 202
	%	72	72			± 25	89	89			± 22

Ref: TEMP-5502 Rev.06