

**LOT**

**BFI 2445**

Rev 1

**CONTROL**



**(Exp.) 2024-03-05**  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS    | BFTROL           |                  |          |          |                         |                  |                  |          |          |                         |
|--------------------------|--------------------|------------------|------------------|----------|----------|-------------------------|------------------|------------------|----------|----------|-------------------------|
|                          |                    | CONTROL          |                  |          |          | TOLERANCES<br>TOLERANCE | CONTROL          |                  |          |          | TOLERANCES<br>TOLERANCE |
|                          |                    | RANGE            |                  | RANGE    |          |                         | RANGE            |                  | RANGE    |          |                         |
|                          |                    | YUMIZEN<br>H1500 | YUMIZEN<br>H2500 | HELO 2.0 | HELO 2.0 |                         | YUMIZEN<br>H1500 | YUMIZEN<br>H2500 | HELO 2.0 | HELO 2.0 |                         |
| BFWBC                    | 10 <sup>6</sup> /l | <b>227</b>       | <b>227</b>       |          |          | ± 91                    | <b>871</b>       | <b>871</b>       |          |          | ± 218                   |
| BFRBC                    | 10 <sup>6</sup> /l | <b>76883</b>     | <b>76883</b>     |          |          | ± 19221                 | <b>498482</b>    | <b>498482</b>    |          |          | ± 124621                |
| BFMN                     | 10 <sup>6</sup> /l | <b>59</b>        | <b>59</b>        |          |          | ± 21                    | <b>87</b>        | <b>87</b>        |          |          | ± 35                    |
|                          | %                  | <b>26</b>        | <b>26</b>        |          |          | ± 9                     | <b>10</b>        | <b>10</b>        |          |          | ± 4                     |
| BFPN                     | 10 <sup>6</sup> /l | <b>168</b>       | <b>168</b>       |          |          | ± 59                    | <b>784</b>       | <b>784</b>       |          |          | ± 196                   |
|                          | %                  | <b>74</b>        | <b>74</b>        |          |          | ± 26                    | <b>90</b>        | <b>90</b>        |          |          | ± 23                    |

Ref: TEMP-5502 Rev.07