

LOT BFI 0123
Rev 1

CONTROL

 (Exp.) 2023-03-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	BFTROL									
		CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE
		RANGE	RANGE	2	RANGE		RANGE	3			
YUMIZEN H1500	YUMIZEN H2500			YUMIZEN H1500	YUMIZEN H2500						
BFWBC	10 ⁶ /l	180	180			± 100	890	890			± 120
BFRBC	10 ⁶ /l	65620	65620			± 12000	459960	459960			± 50000

Ref: TEMP-5502 Rev.05 1300060764