

# Yumizen G DDi 2 for YUMIZEN G1500/G1550



Reference No.:		1300036391	
Lot No.:		931019 (rev. 1)	
Expiry date:		2024-10-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0.2860
		3,00	0.2185
		2,00	0.1350
		1,00	0.0585
		0,50	0.0280
		0,25	0.0140
Slope value		0,07	

# Yumizen G DDi 2 for YUMIZEN G800



Reference No.:		1300036391	
Lot No.:		931019 (rev. 1)	
Expiry date:		2024-10-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0.2665
		3,00	0.2145
		2,00	0.1339
		1,00	0.0535
		0,50	0.0249
		0,25	0.0128
Slope value		0,06	



# Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		931019 (rev. 1)	
Expiry date:		2024-10-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0.1046
		3,00	0.0783
		2,00	0.0540
		1,00	0.0250
		0,50	0.0117
		0,25	0.0062
		Slope value	0,02
<b>Important notice</b>			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type <b>Lin-Lin p-p</b> is selected in <b>Calibration &gt; D-Dimer &gt; ug FEU/ml</b>.          For more information, please contact your local HORIBA Medical representative.</p>			