

# Yumizen G DDi 2 for YUMIZEN G1500/G1550



Reference No.:		1300036391	
Lot No.:		931101 (rev. 1)	
Expiry date:		2024-11-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0.2850
		3,00	0.2170
		2,00	0.1340
		1,00	0.0580
		0,50	0.0275
		0,25	0.0140
Slope value		0,06	

# Yumizen G DDi 2 for YUMIZEN G800



Reference No.:		1300036391	
Lot No.:		931101 (rev. 1)	
Expiry date:		2024-11-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0.2585
		3,00	0.2015
		2,00	0.1280
		1,00	0.0560
		0,50	0.0260
		0,25	0.0130
Slope value		0,06	



# Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		931101 (rev. 1)	
Expiry date:		2024-11-30	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0.0944
		3,00	0.0702
		2,00	0.0508
		1,00	0.0262
		0,50	0.0128
		0,25	0.0065
		Slope value	0,02
<b>Important notice</b>			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type <b>Lin-Lin p-p</b> is selected in <b>Calibration &gt; D-Dimer &gt; ug FEU/ml</b>.          For more information, please contact your local HORIBA Medical representative.</p>			