

ABX Minocal



(Exp.) 2024-03-05

(YYYY - MM - DD)

LOT		CX 491		CAL		(Exp.) 2024-03-05 (YYYY - MM - DD)					
Rev 1		ABX Lyse									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52				4.52					± 0.06
HB HGB	g/dl	13.3				13.3					± 0.2
	g/l	133				133					± 2
	mmol/l	8.26				8.26					± 0.12
HT HCT	%	37.7				38.0					± 1.0
	l/l	0.377				0.380					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	238				239					± 10
VMP MPV	µm ³ ; fl	8.8				9.7					± 0.5

Ref: TEMP-0387 Rev.51 1300090308 FRONT / RECTO