

ABX Minocal

LOT

CX 491

CAL



(Exp.) 2024-03-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8				8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52				4.54					± 0.06
HB HGB	g/dl	13.3				13.4					± 0.2
	g/l	133				134					± 2
	mmol/l	8.26				8.32					± 0.12
HT HCT	%	37.7				38.6					± 1.0
	l/l	0.377				0.386					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	221				241					± 10
VMP MPV	μm ³ ; fl	8.9				9.3					± 0.5

ABX Minocal

LOT

CX 491

Rev 1

CAL



(Exp.) **2024-03-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		Scil Vet ABC Plus ESV 60							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l									± 0.06
HB HGB	g/dl									± 0.2
	g/l									± 2
	mmol/l									± 0.12
HT HCT	%									± 1.0
	l/l									± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l									± 10
VMP MPV	µm ³ ; fl									± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300113555