

ABX Minocal



LOT		CX 491				CAL		 (Exp.) 2024-03-05		(YYYY - MM - DD)		
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0	9.1	9.1	9.2	9.0		9.1	9.1	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.52	4.53	4.52	4.55	4.50		4.52	4.52	± 0.06	
HB HGB	g/dl	13.1	13.1	13.3	13.3	13.2	13.2		13.3	13.3	± 0.2	
	g/l	131	131	133	133	132	132		133	133	± 2	
	mmol/l	8.14	8.14	8.26	8.26	8.20	8.20		8.26	8.26	± 0.12	
HT HCT	%	37.4	37.3	36.5	36.4	36.4	36.0		36.6	36.4	± 1.0	
	l/l	0.374	0.373	0.365	0.364	0.364	0.360		0.366	0.364	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	238	239	227	233	242	236		235	233	± 10	
VMP MPV	µm ³ ; fl	9.0	8.9	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref. TEMP-0387 Rev.51