

ABX Minocal



LOT

CX 491

CAL



(Exp.) 2024-03-05
(YYYY - MM - DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0	9.1	9.1	9.2	9.0		9.1	9.1	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.52	4.53	4.52	4.55	4.50		4.52	4.52	± 0.06
HB HGB	g/dl	13.1	13.1	13.3	13.3	13.2	13.2		13.3	13.3	± 0.2
	g/l	131	131	133	133	132	132		133	133	± 2
	mmol/l	8.14	8.14	8.26	8.26	8.20	8.20		8.26	8.26	± 0.12
HT HCT	%	37.4	37.3	36.5	36.4	36.4	36.0		36.6	36.4	± 1.0
	l/l	0.374	0.373	0.365	0.364	0.364	0.360		0.366	0.364	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	238	239	227	233	242	236		235	233	± 10
VMP MPV	µm ³ ; fl	9.0	8.9	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.51 BACK / VERSO 1300090308