

ABX Minocal



LOT

CX 491

CAL



(Exp.) 2024-03-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75	8.75			8.70	8.70			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57			4.57	4.57			± 0.06
HB HGB	g/dl	12.8	12.8			12.9	12.9			± 0.2
	g/l	128	128			129	129			± 2
	mmol/l	7.95	7.95			8.01	8.01			± 0.12
HT HCT	%	37.5	37.5			36.7	36.3			± 1.0
	l/l	0.375	0.375			0.367	0.363			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	248	248			261	268			± 10
VMP MPV	μm ³ ; fl	10.6	10.6			N/A	N/A			± 0.5

Ref: TEMP-0387 Rev.51 BACK / VERSO 1300113555