

**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT**

**MX 445 L**  
Rev 1

**CONTROL L**

(Exp.)

**2024-04-05**  
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS ES60 Care ST	MICROS	MICROS	MICROS	MICROS 60 8P - 60 18P ADVIA 60	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>2.1</b>				<b>2.1</b>					± 0.4
GR RBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	<b>2.18</b>				<b>2.18</b>					± 0.17
HB HGB	g/dl	<b>5.9</b>				<b>5.9</b>					± 0.5
	g/l	<b>59</b>				<b>59</b>					± 5
	mmol/l	<b>3.66</b>				<b>3.66</b>					± 0.31
HT HCT	%	<b>15.5</b>				<b>15.5</b>					± 2.0
	l/l	<b>0.155</b>				<b>0.155</b>					± 0.020
VGM MCV	µm <sup>3</sup> ; fl	<b>71</b>				<b>71</b>					± 4
TGMH MCH	pg	<b>27.1</b>				<b>27.1</b>					± 2.0
	fmol	<b>1.68</b>				<b>1.68</b>					± 0.12
CCMH MCHC	g/dl	<b>38.1</b>				<b>38.1</b>					± 3.0
	g/l	<b>381</b>				<b>381</b>					± 30
	mmol/l	<b>23.67</b>				<b>23.67</b>					± 1.86
IDR RDW	%	<b>16.0</b>				<b>14.5</b>					± 3.0
PLAQ. PLTS	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>67</b>				<b>65</b>					± 22
VPM MPV	µm <sup>3</sup> ; fl	<b>8.2</b>				<b>7.5</b>					± 2.0
LYMPHO	#	<b>1.39</b>				<b>1.41</b>					± 0.40
	%	<b>66.0</b>				<b>67.0</b>					± 8.0
MONO	#	<b>0.23</b>				<b>0.20</b>					± 0.20
	%	<b>11.0</b>				<b>9.5</b>					± 6.0
GRANULO	#	<b>0.48</b>				<b>0.49</b>					± 0.40
	%	<b>23.0</b>				<b>23.5</b>					± 7.0

**HORIBA ABX SAS**  
 HORIBA ABX SAS  
 Parc Euromédecine  
 Rue du Caducée - BP 7290  
 34184 MONTPELLIER Cedex 4  
 FRANCE  
 Tél.: 33 (0) 4 67 14 15 16  
 Fax.: 33 (0) 4 67 14 15 17

**LOT**

**MX 445 N**  
 Rev 1

**CONTROL N**

(Exp.)

**2024-04-05**  
 (YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST				60 8P - 60 18P ADVIA 60					
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>7.7</b>				<b>7.8</b>					± 0.8
GR RBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	<b>4.52</b>				<b>4.56</b>					± 0.26
HB HGB	g/dl	<b>13.4</b>				<b>13.5</b>					± 0.7
	g/l	<b>134</b>				<b>135</b>					± 7
	mmol/l	<b>8.32</b>				<b>8.38</b>					± 0.43
HT HCT	%	<b>37.1</b>				<b>37.4</b>					± 2.7
	l/l	<b>0.371</b>				<b>0.374</b>					± 0.027
VGM MCV	µm <sup>3</sup> ; fl	<b>82</b>				<b>82</b>					± 4
TGMH MCH	pg	<b>29.6</b>				<b>29.6</b>					± 2.0
	fmol	<b>1.84</b>				<b>1.84</b>					± 0.12
CCMH MCHC	g/dl	<b>36.2</b>				<b>36.1</b>					± 3.0
	g/l	<b>362</b>				<b>361</b>					± 30
	mmol/l	<b>22.45</b>				<b>22.42</b>					± 1.86
IDR RDW	%	<b>14.5</b>				<b>13.5</b>					± 3.0
PLAQ. PLTS	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>260</b>				<b>260</b>					± 40
VPM MPV	µm <sup>3</sup> ; fl	<b>7.5</b>				<b>7.1</b>					± 2.0
LYMPHO	#	<b>2.81</b>				<b>2.89</b>					± 0.40
	%	<b>36.5</b>				<b>37.0</b>					± 5.0
MONO	#	<b>0.50</b>				<b>0.47</b>					± 0.30
	%	<b>6.5</b>				<b>6.0</b>					± 4.0
GRANULO	#	<b>4.39</b>				<b>4.45</b>					± 0.60
	%	<b>57.0</b>				<b>57.0</b>					± 6.0

Ref: TEMP-0830 Rev.42 FRONT / RECTO 1300104144

