

LOT		MX 446 L		Rev 3		CONTROL		L		(Exp.)		2024-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST				60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	2.0				2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.26				2.23						± 0.17	
HB HGB	g/dl	6.1				6.1						± 0.5	
	g/l	61				61						± 5	
	mmol/l	3.79				3.79						± 0.31	
HT HCT	%	16.5				16.3						± 2.0	
	l/l	0.165				0.163						± 0.020	
VGM MCV	µm ³ ; fl	73				73						± 4	
TGMH MCH	pg	27.0				27.4						± 2.0	
	fmol	1.68				1.70						± 0.12	
CCMH MCHC	g/dl	37.0				37.5						± 3.0	
	g/l	370				375						± 30	
	mmol/l	22.96				23.27						± 1.86	
IDR RDW	%	15.0				14.0						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	75				73						± 22	
VPM MPV	µm ³ ; fl	8.6				7.5						± 2.0	
LYMPHO	#	1.27				1.23						± 0.40	
	%	63.3				61.5						± 8.0	
MONO	#	0.20				0.23						± 0.20	
	%	9.6				11.7						± 6.0	
GRANULO	#	0.54				0.54						± 0.40	
	%	27.1				26.8						± 7.0	

LOT		MX 446 N		Rev 3		CONTROL N		(Exp.)		2024-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST				60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.7				7.7					± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.54				4.55					± 0.26
HB HGB	g/dl	13.6				13.6					± 0.7
	g/l	136				136					± 7
	mmol/l	8.45				8.45					± 0.43
HT HCT	%	37.7				37.8					± 2.7
	l/l	0.377				0.378					± 0.027
VGM MCV	µm ³ ; fl	83				83					± 4
TGMH MCH	pg	30.0				29.9					± 2.0
	fmol	1.86				1.86					± 0.12
CCMH MCHC	g/dl	36.1				36.0					± 3.0
	g/l	361				360					± 30
	mmol/l	22.41				22.36					± 1.86
IDR RDW	%	14.5				13.5					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	275				260					± 40
VPM MPV	µm ³ ; fl	7.7				7.5					± 2.0
LYMPHO	#	2.49				2.49					± 0.40
	%	32.4				32.4					± 5.0
MONO	#	1.18				1.18					± 1.06
	%	15.3				15.3					± 13.8
GRANULO	#	4.03				4.03					± 0.77
	%	52.3				52.3					± 10.0

Ref: TEMP-0830 Rev.42 FRONT / RECTO 1300104144

LOT		MX 446 H		Rev 3		CONTROL		H		(Exp.)		2024-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST				60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.5				20.5						± 1.7	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.65				5.68						± 0.32	
HB HGB	g/dl	18.1				18.1						± 0.9	
	g/l	181				181						± 9	
	mmol/l	11.24				11.24						± 0.56	
HT HCT	%	50.9				51.1						± 3.4	
	l/l	0.509				0.511						± 0.034	
VGM MCV	µm ³ ; fl	90				90						± 4	
TGMH MCH	pg	32.0				31.9						± 2.0	
	fmol	1.99				1.98						± 0.12	
CCMH MCHC	g/dl	35.6				35.4						± 3.0	
	g/l	356				354						± 30	
	mmol/l	22.10				21.99						± 1.86	
IDR RDW	%	14.5				13.0						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	550				530						± 65	
VPM MPV	µm ³ ; fl	7.7				7.5						± 2.0	
LYMPHO	#	3.28				3.28						± 1.00	
	%	16.0				16.0						± 5.0	
MONO	#	1.03				1.03						± 0.80	
	%	5.0				5.0						± 4.0	
GRANULO	#	16.20				16.20						± 1.40	
	%	79.0				79.0						± 7.0	