

ABX Minocal



LOT		CX 492		CAL		(Exp.) 2024-04-05 (YYYY - MM - DD)						
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1	9.3	9.2	9.4	9.1		9.2	9.2	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.46	4.48	4.50	4.51	4.48		4.50	4.50	± 0.06	
HB HGB	g/dl	13.2	13.2	13.4	13.4	13.3	13.3		13.4	13.4	± 0.2	
	g/l	132	132	134	134	133	133		134	134	± 2	
	mmol/l	8.20	8.20	8.32	8.32	8.26	8.26		8.32	8.32	± 0.12	
HT HCT	%	37.5	37.2	36.5	36.7	36.5	36.3		36.7	36.7	± 1.0	
	l/l	0.375	0.372	0.365	0.367	0.365	0.363		0.367	0.367	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	239	222	226	235	227		229	226	± 10	
VMP MPV	µm ³ ; fl	8.9	8.9	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref. TEMP-0387 Rev.51