

ABX Minotrol Retic



LOT RX 446
Rev 1

CONTROL

(Exp.) 2024-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500 HELO 1.x				H2500 HELO 1.x				H2500 HELO 1.x			
Méthode Manuelle / Manual													
RET	%	1.6			± 0.8	5.4			± 2.2	10.4			± 3.2
RET	%	1.63			± 0.82	4.71			± 1.18	9.41			± 1.88
RET	10 ⁶ /mm ³	0.057			± 0.029	0.165			± 0.041	0.323			± 0.065
	10 ⁹ /l	57			± 29	165			± 41	323			± 65
MFI MRI *	%	23.4			± 8.0	27.4			± 10.0	28.5			± 12.0
	-	0.234			± 0.080	0.274			± 0.100	0.285			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.51			± 0.30	3.50			± 0.30	3.43			± 0.30
PIC RANGE	-	19			± 8	19			± 8	19			± 8

Ref: TEMP-0382 Rev.45 FRONT / RECTO 1300116193

* For research use only in Canada and USA