

ABX **Micros ES** 60

ABX **Micros ESV** 60

**Micros Care ST**

**scil Vet abc Plus<sup>+</sup>**

## Printers List

Ref.: RAL095F

## Printers Validated by HORIBA Medical

The compatibility of the printers listed below has been tested. HORIBA Medical can guarantee the functioning of ABX Micros ES60 OT/CT, ABX Micros ESV60, Micros Care ST and scil Vet abc Plus+ instruments with these printers:

Printer	HORIBA Medical part number for ordering
EPSON LX 350 110V (USB)	1300011096
EPSON LX 350 220V (USB)	1300011095
EPSON AL-M320DN 230V	1300079681
EPSON M300 220V	1209101727

## Other Printers

Other printers (not tested by HORIBA Medical) can be used with the drivers integrated in the instruments.

To identify the printers and the compatible drivers, visit the Linux Foundation website:

<https://openprinting.org/printers>.



Some printers may require purchasing a USB cord.

HORIBA Medical recommends using printers compatible with protocol IPP (Internet Printing Protocol).

---