



Manufactured for  
 HORIBA Instruments Incorporated:  
 Pointe Brand  
 5449 Research Drive  
 Canton, MI 48188



**LOT** 319307

**SAP** 5390015524



**REF** T7535-CAL



2024-09-30

**Calibrator Values\***

	Valeur titrée Assigned value Deklarierter Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tildelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet
<b>Calibrator A</b>	0	mg/dL
	0.00	g/L
<b>Calibrator B</b>	77	mg/dL
	0.77	g/L
<b>Calibrator C</b>	155	mg/dL
	1.55	g/L
<b>Calibrator D</b>	258	mg/dL
	2.58	g/L
<b>Calibrator E</b>	361	mg/dL
	3.61	g/L
<b>Calibrator F</b>	516	mg/dL
	5.16	g/L

\*Standardized against CAP/BCR/IFCC Reference Preparation for Proteins in Human Serum (RPPHS), lot 91/0619.