



Manufactured for
HORIBA Instruments Incorporated:
Pointe Brand
5449 Research Drive
Canton, MI 48188



LOT 319303

REF F7505-CTL

SAP 5390015523



2024-05-31

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet	Level I Exp. Date: 2024-05-31 Valeur titrée Assigned value Deklarierter Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tilldelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd			Level II Exp. Date: 2024-05-31 Valeur titrée Assigned value Deklarierter Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tilldelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd		
		Mean	SD	Range	Mean	SD	Range
Ferritin (1st I.S. Standardization)	ng/mL	101	7.5	86-116	396	29.5	337-455
Ferritin (3rd I.S. Standardization)	ng/mL	120	9.0	102-138	471	35.5	400-542

Ferritin values are provided with two options: 1st or 3rd I.S. (International Standard) Values.

It is recommended to use the same standardized system as was chosen when selecting values for the Pointe Ferritin Calibrators (F7505-CAL).