

LOT BFI 2446
Rev 1

CONTROL

 (Exp.) 2024-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	BFTROL										
		CONTROL					TOLERANCES TOLERANCE	CONTROL				
		RANGE		RANGE		RANGE		RANGE		TOLERANCES TOLERANCE		
		YUMIZEN H1500	YUMIZEN H2500	HELO 2.0	HELO 2.0	YUMIZEN H1500		YUMIZEN H2500	HELO 2.0		HELO 2.0	
BFWBC	10 ⁶ /l	231	231			± 92	852	852			± 213	
BFRBC	10 ⁶ /l	76800	76800			± 19200	503561	503561			± 125890	
BFMN	10 ⁶ /l	69	69			± 24	108	108			± 43	
	%	30	30			± 11	13	13			± 5	
BFPN	10 ⁶ /l	163	163			± 57	744	744			± 186	
	%	70	70			± 25	87	87			± 22	

Ref: TEMP-5502 Rev.07