

ABX Minocal



(Exp.) 2024-05-05

(YYYY - MM - DD)

LOT		CX 493		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST				60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1				9.0						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48				4.48						± 0.06
HB HGB	g/dl	13.4				13.5						± 0.2
	g/l	134				135						± 2
	mmol/l	8.32				8.38						± 0.12
HT HCT	%	38.1				38.1						± 1.0
	l/l	0.381				0.381						± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	236				232						± 10
VMP MPV	μm ³ ; fl	8.7				7.6						± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300090308