

ABX Minocal



LOT

CX 493

CAL



(Exp.) 2024-05-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1				9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48				4.48					± 0.06
HB HGB	g/dl	13.4				13.5					± 0.2
	g/l	134				135					± 2
	mmol/l	8.32				8.38					± 0.12
HT HCT	%	38.1				38.1					± 1.0
	l/l	0.381				0.381					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	236				232					± 10
VMP MPV	µm ³ ; fl	8.7				7.6					± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300090308