

ABX Minocal



LOT		CX 493		CAL		(Exp.) 2024-05-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48				4.43					± 0.06
HB HGB	g/dl	13.6				13.4					± 0.2
	g/l	136				134					± 2
	mmol/l	8.45				8.32					± 0.12
HT HCT	%	38.1				38.5					± 1.0
	l/l	0.381				0.385					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	224				235					± 10
VMP MPV	µm ³ ; fl	9.2				9.2					± 0.5

ABX Minocal



LOT		CX 493		CAL		 (Exp.) 2024-05-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
			Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.8								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.46								± 0.06
HB HGB	g/dl		13.3								± 0.2
	g/l		133								± 2
	mmol/l		8.26								± 0.12
HT HCT	%		37.9								± 1.0
	l/l		0.379								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		225								± 10
VMP MPV	µm ³ ; fl		8.6								± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300113555