

**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT**

**MC 431**

**CONTROL**



(Exp.) **2024-07-05**  
(YYYY - MM - DD)

Rev 2

		ABX Lyse															
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				1	CONTROL				2	CONTROL				3	TOLERANCES TOLERANCE
		MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	
				60 8P 60 18P					60 8P 60 18P					60 8P 60 18P			
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l			<b>2.1</b>		± 0.4			<b>7.9</b>		± 1.0			<b>21.1</b>		± 2.2	
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l			<b>2.25</b>		± 0.15			<b>4.54</b>		± 0.18			<b>5.59</b>		± 0.25	
HB HGB	g/dl			<b>6.2</b>		± 0.4			<b>13.4</b>		± 0.5			<b>17.7</b>		± 0.7	
	g/l			<b>62</b>		± 4			<b>134</b>		± 5			<b>177</b>		± 7	
	mmol/l			<b>3.85</b>		± 0.25			<b>8.32</b>		± 0.31			<b>10.99</b>		± 0.44	
HT HCT	%			<b>16.7</b>		± 2.0			<b>36.8</b>		± 2.5			<b>49.8</b>		± 3.0	
	l/l			<b>0.167</b>		± 0.020			<b>0.368</b>		± 0.025			<b>0.498</b>		± 0.030	
VGM MCV	µm <sup>3</sup> ; fl			<b>74</b>		± 4			<b>81</b>		± 4			<b>89</b>		± 4	
TGMH MCH	pg			<b>27.6</b>		± 2.4			<b>29.5</b>		± 2.4			<b>31.7</b>		± 2.8	
	fmol			<b>1.71</b>		± 0.15			<b>1.83</b>		± 0.15			<b>1.97</b>		± 0.18	
CCMH MCHC	g/dl			<b>37.2</b>		± 3.0			<b>36.4</b>		± 3.0			<b>35.6</b>		± 3.0	
	g/l			<b>372</b>		± 30			<b>364</b>		± 30			<b>356</b>		± 30	
	mmol/l			<b>23.12</b>		± 1.86			<b>22.63</b>		± 1.86			<b>22.09</b>		± 1.86	
IDR RDW	%			<b>14.0</b>		± 3.0			<b>14.0</b>		± 3.0			<b>13.0</b>		± 3.0	
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l			<b>80</b>		± 20			<b>268</b>		± 40			<b>532</b>		± 60	
VPM MPV	µm <sup>3</sup> ; fl			<b>7.5</b>		± 2.0			<b>6.8</b>		± 2.0			<b>7.1</b>		± 2.0	
LYMPHO	#			<b>1.3</b>		± 0.2			<b>2.7</b>		± 0.5			<b>3.3</b>		± 1.1	
	%			<b>64.0</b>		± 8.0			<b>34.0</b>		± 6.0			<b>15.5</b>		± 5.0	
MONO	#			<b>0.2</b>		± 0.2			<b>1.2</b>		± 1.1			<b>1.7</b>		± 1.3	
	%			<b>10.0</b>		± 7.0			<b>15.3</b>		± 13.8			<b>8.0</b>		± 6.0	
GRANULO	#			<b>0.5</b>		± 0.2			<b>4.0</b>		± 1.0			<b>16.1</b>		± 2.1	
	%			<b>26.0</b>		± 7.0			<b>50.7</b>		± 10.0			<b>76.5</b>		± 10.0	
CRP	mg/dl			N/A		± 0.30			N/A		± 0.70			N/A		± 1.00	
	mg/l			N/A		± 3.0			N/A		± 7.0			N/A		± 10.0	

Ref: TEMP-1070 Rev.42 FRONT/RECTO 1300092170