

ABX Minocal

LOT

CX 493

CAL



(Exp.) 2024-05-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48				4.43					± 0.06
HB HGB	g/dl	13.6				13.4					± 0.2
	g/l	136				134					± 2
	mmol/l	8.45				8.32					± 0.12
HT HCT	%	38.1				38.5					± 1.0
	l/l	0.381				0.385					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	224				235					± 10
VMP MPV	µm ³ ; fl	9.2				9.2					± 0.5

ABX Minocal

LOT

CX 493

Rev 1

CAL



(Exp.) **2024-05-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		Scil Vet ABC Plus ESV 60							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l									± 0.06
HB HGB	g/dl									± 0.2
	g/l									± 2
	mmol/l									± 0.12
HT HCT	%									± 1.0
	l/l									± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l									± 10
VMP MPV	µm ³ ; fl									± 0.5

Ref: TEMP-0387 Rev.51 FRONT / RECTO 1300113555