

ABX Minocal



LOT		CX 493		CAL		 (Exp.) 2024-05-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1	9.2	9.0	9.3	9.0		9.1	9.0	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.44	4.44	4.41	4.46	4.48	4.42		4.47	4.46	± 0.06	
HB HGB	g/dl	13.3	13.2	13.3	13.5	13.4	13.4		13.5	13.5	± 0.2	
	g/l	133	132	133	135	134	134		135	135	± 2	
	mmol/l	8.26	8.20	8.26	8.38	8.32	8.32		8.38	8.38	± 0.12	
HT HCT	%	37.5	37.5	36.2	37.0	36.5	36.2		37.1	37.0	± 1.0	
	l/l	0.375	0.375	0.362	0.370	0.365	0.362		0.371	0.370	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	232	234	227	228	235	230		231	228	± 10	
VMP MPV	µm ³ ; fl	8.8	9.0	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref. TEMP-0387 Rev.51