

# ABX Minocal



| LOT                      |                                                       | CX 493  |             | CAL    |        | <div><div></div><div>(Exp.) 2024-05-05</div><div>(YYYY - MM - DD)</div></div> |        |        |                     |        |                         |  |
|--------------------------|-------------------------------------------------------|---------|-------------|--------|--------|-------------------------------------------------------------------------------|--------|--------|---------------------|--------|-------------------------|--|
|                          |                                                       | Rev 1   | ABX Lysebio |        |        |                                                                               |        |        |                     |        |                         |  |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | MICROS  | MICROS      | PENTRA | PENTRA | PENTRA                                                                        | PENTRA | PENTRA | PENTRA              | PENTRA | TOLERANCES<br>TOLERANCE |  |
|                          |                                                       | EMI CRP | LC-767G     | MS CRP | XLR    | DX NEXUS<br>DF NEXUS                                                          | MS60   |        | 60<br>60 C+<br>ES60 | XL80   |                         |  |
|                          |                                                       |         |             |        |        |                                                                               |        |        |                     |        |                         |  |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.1     | 9.1         | 9.2    | 9.0    | 9.3                                                                           | 9.0    |        | 9.1                 | 9.0    | ± 0.2                   |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.44    | 4.44        | 4.41   | 4.46   | 4.48                                                                          | 4.42   |        | 4.47                | 4.46   | ± 0.06                  |  |
| HB HGB                   | g/dl                                                  | 13.3    | 13.2        | 13.3   | 13.5   | 13.4                                                                          | 13.4   |        | 13.5                | 13.5   | ± 0.2                   |  |
|                          | g/l                                                   | 133     | 132         | 133    | 135    | 134                                                                           | 134    |        | 135                 | 135    | ± 2                     |  |
|                          | mmol/l                                                | 8.26    | 8.20        | 8.26   | 8.38   | 8.32                                                                          | 8.32   |        | 8.38                | 8.38   | ± 0.12                  |  |
| HT HCT                   | %                                                     | 37.5    | 37.5        | 36.2   | 37.0   | 36.5                                                                          | 36.2   |        | 37.1                | 37.0   | ± 1.0                   |  |
|                          | l/l                                                   | 0.375   | 0.375       | 0.362  | 0.370  | 0.365                                                                         | 0.362  |        | 0.371               | 0.370  | ± 0.010                 |  |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 232     | 234         | 227    | 228    | 235                                                                           | 230    |        | 231                 | 228    | ± 10                    |  |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.8     | 9.0         | N/A    | N/A    | N/A                                                                           | N/A    |        | N/A                 | N/A    | ± 0.5                   |  |

1300090308  
BACK / VERSO  
Ref. TEMP-0387 Rev.51