

# ABX Minocal



**LOT**

**CX 493**

**CAL**



**(Exp.) 2024-05-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500	H1500	H2500	H2500					
		HELO 1.x	HELO 2.0							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.65</b>	<b>8.65</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.53</b>	<b>4.56</b>							± 0.06
HB HGB	g/dl	<b>13.5</b>	<b>13.5</b>							± 0.2
	g/l	<b>135</b>	<b>135</b>							± 2
	mmol/l	<b>8.38</b>	<b>8.38</b>							± 0.12
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>240</b>	<b>238</b>							± 10
VGM MCV	μm <sup>3</sup> ; fl	<b>83.0</b>	<b>83.9</b>							± 4.0

1300116192  
Ref: TEMP-0387 Rev.51