

ABX Minocal



LOT

CX 493

CAL



(Exp.) 2024-05-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.95	8.95			8.75	8.75			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51			4.51	4.51			± 0.06
HB HGB	g/dl	13.1	13.1			13.2	13.2			± 0.2
	g/l	131	131			132	132			± 2
	mmol/l	8.14	8.14			8.20	8.20			± 0.12
HT HCT	%	37.9	37.9			36.4	36.6			± 1.0
	l/l	0.379	0.379			0.364	0.366			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	236	236			239	260			± 10
VMP MPV	μm ³ ; fl	10.1	10.1			N/A	N/A			± 0.5

Ref: TEMP-0387 Rev.51 BACK / VERSO 1300113555