

LOT		MX 447 L		CONTROL		L		(Exp.)		2024-08-05 (YYYY-MM-DD)	
Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS								TOLERANCES TOLERANCE	
		Scil Vet ABC Plus Scil Vet ABC ESV 60									
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.8									± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.27									± 0.15
HB HGB	g/dl	5.9									± 0.4
	g/l	59									± 4
	mmol/l	3.66									± 0.25
HT HCT	%	16.3									± 2.0
	l/l	0.163									± 0.020
VGM MCV	µm ³ ; fl	72									± 4
TGMH MCH	pg	26.0									± 2.0
	fmol	1.61									± 0.12
CCMH MCHC	g/dl	36.1									± 3.0
	g/l	361									± 30
	mmol/l	22.42									± 1.86
IDR RDW *	%	15.0									± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	74									± 20
VPM MPV *	µm ³ ; fl	8.0									± 2.0
LYMPHO *	#	0.96									± 0.40
	%	53.5									± 8.0
MONO *	#	0.23									± 0.20
	%	13.0									± 6.0
GRANULO *	#	0.60									± 0.40
	%	33.5									± 7.0

* N/A for Micros Scil Vet ABC

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT		MX 447 N	CONTROL		N					(Exp.)	2024-08-05 (YYYY-MM-DD)	
		Rev 1	LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS									TOLERANCES TOLERANCE	
			Scil Vet ABC Plus Scil Vet ABC ESV 60									
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.1								± 0.8	
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l		4.54								± 0.18	
HB HGB	g/dl		13.2								± 0.5	
	g/l		132								± 5	
	mmol/l		8.20								± 0.31	
HT HCT	%		36.8								± 2.5	
	l/l		0.368								± 0.025	
VGM MCV	µm ³ ; fl		81								± 4	
TGMH MCH	pg		29.1								± 2.0	
	fmol		1.81								± 0.12	
CCMH MCHC	g/dl		35.9								± 3.0	
	g/l		359								± 30	
	mmol/l		22.29								± 1.86	
IDR RDW *	%		14.5								± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		262								± 40	
VPM MPV *	µm ³ ; fl		7.5								± 2.0	
LYMPHO *	#		2.02								± 0.40	
	%		28.5								± 5.0	
MONO *	#		0.71								± 0.45	
	%		10.0								± 6.0	
GRANULO *	#		4.37								± 0.60	
	%		61.5								± 6.0	

* N/A for Micros Scil Vet ABC

1300104143 BACK / VERSO Ref: TEMP-0630 Rev.43

LOT

MX 447 H
Rev 1

CONTROL

H

 (Exp.)

2024-08-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS								
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		19.7							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l		5.73							± 0.20
HB HGB	g/dl		18.1							± 0.6
	g/l		181							± 6
	mmol/l		11.24							± 0.37
HT HCT	%		50.4							± 3.0
	l/l		0.504							± 0.030
VGM MCV	µm ³ ; fl		88							± 4
TGMH MCH	pg		31.6							± 2.0
	fmol		1.96							± 0.12
CCMH MCHC	g/dl		35.9							± 3.0
	g/l		359							± 30
	mmol/l		22.29							± 1.86
IDR RDW *	%		13.5							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		498							± 45
VPM MPV *	µm ³ ; fl		7.4							± 2.0
LYMPHO *	#		2.17							± 1.00
	%		11.0							± 5.0
MONO *	#		1.38							± 0.80
	%		7.0							± 4.0
GRANULO *	#		16.15							± 1.40
	%		82.0							± 7.0

* N/A for Micros Scil Vet ABC