

ABX Erytrol



LOT

NX 447

Rev 1

CONTROL



(Exp.) 2024-07-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
			PENTRA				PENTRA				PENTRA		
	DX NEXUS			DX NEXUS			DX NEXUS			DX NEXUS			
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.5		± 0.8		18.5		± 1.5		2.6		± 0.4
GBC CWBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.2		± 0.8		16.8		± 1.4		1.5		± 0.4
ERB	%		4.2		± 1.7		10.0		± 3.3		78.7		± 50.0
ERB	#		0.302		± 0.100		1.680		± 0.700		1.181		± 0.500

1219930076-C

Ref: TEMP-0822 Rev.37

GBC / CWBC

* les valeurs sont données à titre d'indication et doivent être utilisées comme RuO (Utilisation en Recherche Seulement) Indicative expected values, must be used as RuO (Research used only)