

ABX Minocal



LOT		CX 494		CAL		(Exp.) 2024-06-05 (YYYY - MM - DD)						
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1	9.4	9.1	9.3	9.1		9.2	9.1	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.45	4.35	4.45	4.46	4.43		4.48	4.45	± 0.06	
HB HGB	g/dl	13.1	13.2	13.1	13.3	13.2	13.2		13.2	13.3	± 0.2	
	g/l	131	132	131	133	132	132		132	133	± 2	
	mmol/l	8.14	8.20	8.14	8.26	8.20	8.20		8.20	8.26	± 0.12	
HT HCT	%	36.9	37.2	35.5	36.5	35.9	35.9		36.5	36.5	± 1.0	
	l/l	0.369	0.372	0.355	0.365	0.359	0.359		0.365	0.365	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	239	236	221	238	241	232		240	238	± 10	
VMP MPV	µm ³ ; fl	8.6	8.6	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref. TEMP-0387 Rev.51